



Agilent Technologies

## Certificate of Analysis

### ESI-L Low Concentration Tuning Mix

Agilent Part Number: G1969-85000

Sample Lot Number: CS-3025

#### Concentration (weight to volume) and Purity/Grades:

Neat Material	CAS #	Gravimetric Conc.	Purity
Betaine	107-43-7	<0.01%	99.3%
Trifluoroacetic acid ammonium salt	3336-58-1	<0.01%	100%
Hexamethoxyphosphazine	957-13-1	<0.01%	99.0%
Hexakis(2,2-difluoroethoxy)phosphazine	186817-57-2	<0.01%	99.0%
Hexakis(1H, 1H, 3H-tetrafluoropropoxy)phosphazine	58943-98-9	<0.01%	99.0%
Hexakis(1H, 1H, 5H-octafluoropentoxy)phosphazine	16059-16-8	<0.01%	97.0%
Hexakis(1H, 1H, 7H-dodecafluoroheptoxy)phosphazine	3830-74-8	<0.01%	99.0%
Hexakis(1H, 1H, 9H-perfluorononyloxy)phosphazine	186043-67-4	<0.01%	99.0%
Hexakis(1H, 1H, 4H-hexafluorobutyloxy)phosphazine	186406-47-3	<0.01%	97.0%
Hexakis(1H, 1H, 6H-decafluorohexyloxy)phosphazine	186406-48-4	<0.01%	97.0%
Hexakis(1H, 1H, 8H-tetradecafluorooctyloxy)phosphazine	186406-49-5	<0.01%	98.0%
Tris(trifluoromethyl)-1, 3, 5-triazine	368-66-1	<0.01%	99.8%
Tris(heptafluoropropyl)-1, 3, 5-triazine	915-76-4	<0.01%	99.8%
<b>Solvent:</b> Acetonitrile (HPLC Grade)	75-05-8	95.0%	99.9%
DI Water	7732-18-5	5.0%	De-ionized

**Storage Condition:** Store Refrigerated (2° to 8° C)

#### Traceability:


This standard has been produced gravimetrically using ISO9001 quality procedures. NIST traceable weights are used to verify balance calibration with the preparation of each lot. Concentration of analyte in solution is ug/ml +/- 0.5%, uncertainty based upon balance and Class A volumetric glassware. API-Mass spectrometry was used to evaluate this multi-standard solution.



## G1969-85000 POS ES-TOF SPECTRA

☐ +TOF MS: 0.097 to 1.967 min from QC100203\_pp6\_ESpos.wiff Agilent

Max. 4.0e5 counts.



Sample: G1969-85000

Instrument: Agilent G1969A API-TOF MS

Mode:ES

Polarity: Positive ion mode

CDS infusion rate: 0.1 ml/min

Vcap: -4 kV

Fragmentor: 215 V

Nebulizer pressure : 20 psig

Drying gas flow : 6 L/min

Drying gas temperature : 300 °C

MCP : 650 V

PMT : 703 V

Scan range: m/z 50-3200


Transients : 10,000

Flight tube: -6.5 kV

Oct. RF : 250 V pk



## G1969-85000 NEG ES-TOF SPECTRA



Sample: G1969-85000 (1:4 dilution in 95:5 ACN/H<sub>2</sub>O)

Instrument: Agilent G1969A API-TOF MS

Mode:ES

Polarity: Negative ion mode

CDS infusion rate: 0.1 ml/min

Vcap: +4 kV

Fragmentor: 160 V

Nebulizer pressure : 20 psig

Drying gas flow : 6 L/min

Drying gas temperature : 300 °C

MCP : 650 V

PMT : 703 V

Scan range: m/z 50-3200

Transients : 10,000

Flight tube: +6.5 kV

Oct. RF : 250 V pk




**PRINCIPAL IONS**

<b>MASS</b>	<b>POS</b>	<b>NEG</b>
<b>1</b>	<b>118.086255</b>	<b>112.985587</b>
<b>2</b>	<b>322.048121</b>	<b>301.998139</b>
<b>3</b>	<b>622.028960</b>	<b>601.978977</b>
<b>4</b>	<b>922.009798</b>	<b>1033.988109</b>
<b>5</b>	<b>1221.990637</b>	<b>1333.968947</b>
<b>6</b>	<b>1521.971475</b>	<b>1633.949786</b>
<b>7</b>	<b>1821.952313</b>	<b>1933.930624</b>
<b>8</b>	<b>2121.933152</b>	<b>2233.911463</b>
<b>9</b>	<b>2421.913990</b>	<b>2533.892301</b>
<b>10</b>	<b>2721.894829</b>	<b>2833.873139</b>

**Date of Release:** 01 June 2018

**Date of Expiration:** 30 June 2020

  
John Russo  
President

  
Monica Bourgeois  
Director of QA/RA