

Tel: (800)356-1688
Fax: (814)353-1309

www.restek.com

Gravimetric Certificate

FOR LABORATORY USE ONLY-READ SDS PRIOR TO USE.

This Reference Material is intended for Laboratory Use Only as a standard for the qualitative and/or quantitative determination of the analyte(s) listed.

Catalog No. : 31858 Lot No.: A0144705

Description: n-Decane

n-Decane (neat), 1mL/ampul

Container Size: 2 mL Pkg Amt: > 1 mL

Compound #	Compound	CAS#	Percent Purity	Actual Weight	Expanded Uncertainty (95% C.L.; K=2)
1	n-Decane (C10)	124-18-5	99%	1.0 ml	+/- 0.0058 ml
Solvent:	None				

1. i. f. Carat

Cydnei I Crust - Mix Technician

Date Mixed: 07-Jan-2019

Balance: B251644995

Manufactured under Restek's ISO 9001:2008 Registered Quality System Certificate #FM 80397

General Reference Material Notes

Expiration Notes:

- Expiration date valid for unopened ampul stored in compliance with the recommended conditions.
- Uncertainty, concentration, and expiration of the RM are based on the unopened product being stored according to the recommended condition found in the storage field.

Purity Notes:

- Purity and/or chemical identity are determined by one or more of the following techniques: GC/FID, HPLC, GC/µECD, GC/MS, LC/MS, RI, and/or melting point.
- Compounds with a listed purity of less than 99% have been weight corrected to compensate for impurities and/or salts.
 A correction factor is used to calculate the amount of compound necessary to achieve the desired concentration of the parent compound in solution.
- Purity of isomeric compounds is reported as the sum of the isomers.
- Purity values are rounded to the nearest whole number.

Uncertainty Value Notes:

• Uncertainties are determined using data from balances and glassware, raw material purity, and, when significant, equipment tolerances or calibration results.

Manufacturing Notes:

• Concentration is based upon gravimetric preparation using either a balance whose calibration has been verified daily using NIST traceable weights, and/or dilutions with Class A glassware.

Handling Notes:

- Stability of the unopened product, when stored in compliance with the recommended conditions, is guaranteed through
 the expiration displayed on the product label and certificate. Contact Restek for additional opened product stability
 information, with the knowledge/understanding that open product stability is subject to the specific handling and
 environmental conditions to which the product is exposed. For your convenience Restek supplies deactivated vials with
 most standards packed in 2mL ampules. Larger volume deactivated vials are available through Restek as a custom
 ordered item. Additionally, Restek sells DMDCS for the purpose of glassware deactivation as catalog number 31861,
 which includes complete instructions.
- If any undissolved material is visible inside the ampul, sonicate the unopened ampul until the material is completely dissolved.